

# Planmeca Chair™

Excellent  
working ergonomics  
for the dental team  
and optimal comfort  
for patients of  
all sizes.



# Intelligent design details

Ergonomic in its design, the **Planmeca Chair™** offers ultimate comfort for both the patient and the dental team.



## Perfect patient comfort

The **Planmeca Chair™** is available either with a fixed or automatic legrest. The *fixed* legrest provides optimal knee angulation and ensures uniform upholstery. Patient entry and exit can be facilitated with the *automatic* foldable legrest. The perfectly designed backrest and firm armrests maximise patient comfort.

## Easily adjustable headrest

The double articulated headrest allows comfortable positioning of the patient's head and also provides good visibility of the treatment area. Due to the smooth and lightweight headrest length adjustment, the chair can be easily adapted for patients of all sizes.



## Excellent working ergonomics

The slim and narrow backrest and extensive legroom ensure optimal working ergonomics for the whole dental team. The automatic legrest facilitates closer face-to-face positioning that is greatly appreciated in doctor-patient consultations and in prosthodontic treatments.



# Options for a variety of needs

The highly functional **Planmecca Chair™** adapts perfectly to a variety of specialised needs.

## Flexible positioning

The 180-degree swivel allows the chair to be swiveled 90 degrees left and right. This is a valuable feature when using an intraoral X-ray unit or other auxiliary equipment. The swivel function provides a unique freedom for the design of the clinic and is especially appreciated when space is limited. The mobile base with wheels provides extra possibilities for positioning the Planmecca Chair.

## Great adaptability

The wide movement range enables treating patients freely in various positions. The chair can be adjusted to a very low position, a useful feature when the patient is in a sitting or semi-reclined position. The chair can also be adjusted to a high position if the dentist prefers to perform a treatment in a standing position.



## Fit for surgical needs

A surgical intravenous armrest is available for immobilising the patient's arm for example during sedation.

The surgical headrest option provides enhanced sideways support for the patient's head.



## Unbeatable combination for specialist treatments

The **Planmecca Chair™** can be equipped with the unbeatable combination of our large over-the-patient tray and the **Planmecca SingLED®** operating light. This makes the chair an ideal choice for both orthodontic and surgical treatments.

The large and solid over-the-patient tray is made from cast aluminium and has an outstanding load-bearing capacity.

The Planmecca SingLED operating light provides optimised illumination of the treatment area and can be easily adjusted from the upper to the lower jaw position.

# High-quality upholsteries for enduring comfort



The **Comfy™** upholstery is made of durable artificial leather. The thin shape ensures excellent reach to the treatment area and the seamless design allows easy upkeep and enhanced infection control.






Made of extremely durable artificial leather and viscoelastic memory foam, the **Ultra Relax™** upholstery adapts perfectly to patient's body shape, allowing them to stay relaxed even during the longest treatment sessions.

Both **Comfy™** and **Ultra Relax™** upholsteries come in a wide range of stylish colours.

# Technical specifications

Power requirements	100 VAC, 115 VAC, 220–240 VAC
Mains frequency	50/60 Hz
Swivel	±90°, total 180°
Weight	120 kg (290 lbs)
Mounting	To be fixed to the floor Mobile version as an option
Frame	Cast aluminium alloy
Colour	White (RAL 9016)



The fixed legrest has optimum knee angulation.

The optional automatic legrest can be used when the chair is swivelled 90° to the right or left.

# Stand out with colour



## PLANMECA

Asentajankatu 6 | 00880 Helsinki | Finland | tel. +358 20 7795 500 | sales@planmeca.com | www.planmeca.com

Images may contain optional items not included in standard delivery. Available configurations and features may have country or area specific variations. Some products displayed above may not be available in all countries or areas. Rights for changes reserved.

Planmeca, All in one, Anatomat Plus, Cobra, Comfy, DentreVac, Digital perfection, Economat Plus, Elegant, Flexy, Mini-dent, Perio Fresh, PlanEasyMill, Planmeca 4D, Planmeca AINO, Planmeca ARA, Planmeca CAD/CAM, Planmeca CALM, Planmeca Chair, Planmeca Clarify, Planmeca Compact, Planmeca Creo, Planmeca Emerald, Planmeca FIT, Planmeca Intra, Planmeca iRomexis, Planmeca Lumion, Planmeca Lumo, Planmeca Maximity, Planmeca Minea, Planmeca Minendo, Planmeca Minetto, Planmeca mRomexis, Planmeca Noma, Planmeca Olo, Planmeca Online, Planmeca PlanCAD, Planmeca PlanCAM, Planmeca PlanClear, Planmeca PlanID, Planmeca PlanMill, Planmeca Planosil, Planmeca PlanPure, Planmeca PlanScan, Planmeca PlanView, Planmeca ProCeph, Planmeca ProFace, Planmeca ProID, Planmeca ProMax, Planmeca ProModel, Planmeca ProOne, Planmeca ProScanner, Planmeca ProSensor, Planmeca ProX, Planmeca Romexis, Planmeca Serenus, Planmeca SingLED, Planmeca SmartGUI, Planmeca Solea, Planmeca Sovereign, Planmeca Ultra Low Dose, Planmeca Vision, Planmeca Viso, Planmeca Verity, Planmeca Waterline Cleaning System, Planmeca Xtremity, Proline Dental Stool, ProTouch, Saddle Stool, SmartPan, SmartTouch, Trendy, and Ultra Relax are registered or non-registered trademarks of Planmeca in various countries.

